

JRWD503 Power Wheelchair



Features

- Dual Functionality allows the wheelchair to operate as an electric wheelchair or as a manual wheelchair.
- Raising the armrests allows you to get closer to tables and makes it easier to transfer.
- Detachable battery pack and foldable frame to make it easier to transport in a car
- Low-cost, yet powerful and reliable.
- 20AH Lead-acid batteries are low-cost to replace and give 15 – 20 km or travel range.
- Safety Mark charger provided.



Dual Functionality

Operate as an electric wheelchair or as a manual wheelchair



Flip-up Armrest

For easy chair-to-chair sideways transfer



Compact Design

Allows for easy folding and storage in small spaces

Freedom on Wheels

FALCON
MOBILITY

SPECIFICATIONS

JRWD503 DUAL FUNCTION FOLDABLE POWER WHEELCHAIR

ITEM #	JRWD503
Overall width	66.5 cm
Overall length	108 cm
Overall height	94 cm
Frame Weight	30 kg
Battery Weight	14 kg
Overall Weight	44 kg
Seat Width	40 cm or 45 cm (16" or 18")
Seat Depth	43 cm (17")
Backrest Height	38 cm (15")
Seat-to-Floor Height (Back)	51 cm (20")
Seat-to-Floor Height (Back)	49 cm *19")
Footrest Type	Swingaway, Detachable
Footrest Length	14" to 16", adjustable in 1" increments
Armrest Type	Flip-up
Seat-to-Armrest Height	9"
Armrest Size	13" x 2"
Driving Range	15 – 20 km
Max Speed	8 km/h
Turning Radius	75 cm" (30")
Rear Wheel Diameter	24"
Front Wheel Diameter	8"
Maximum Capacity	120 kg
Batteries	2 x 12V 20AH SLA
Charger	2A c/w Safety Mark
Motors	250w
Warranty	One year on electronics, motors and frame 6 months on batteries

DISTRIBUTED BY:



Head Office (☎: +65 6774 7518)

140 PAYA LEBAR ROAD
#05-11 AZ @ PAYA LEBAR
SINGAPORE 409015

OPENING HOURS:
MON – FRI: 10.00 AM to 6.00 PM
SAT : 10.30 AM to 4.30 PM

🌐: www.falconmobility.com.sg

✉: sales@falconmobility.com.sg